






Clinical Pathology Tube Guide

IDEXX BioAnalytics provides these tubes at no charge. To request tubes or for any other questions contact the Research Support Team at 1-800-669-0825 or idxxbioanalytics@idexx.com



Sample Type	Serum* Blood Chemistry, Immunology, Endocrinology	Whole Blood Hematology, Fluids or Cytology
Tube Types	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>SST-SERUM Separator Tube 4mL or 8mL</p> <p>No anticoagulant red cap with yellow ring, clot activator with gel</p> </div> <div style="text-align: center;">  <p>RTT-RED Top Tube 4 mL or 9 mL</p> <p>Black ring, with clot activator</p> </div> <div style="text-align: center;">  <p>SST-MINI Gold Top Tube 0.8 mL</p> <p>Clot activator and gel</p> </div> </div>	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>LTT-LAVENDER Top Tube 2mL or 9mL</p> <p>EDTA anticoagulant lavender cap with white ring</p> </div> <div style="text-align: center;">  <p>LTT-MINI Top Tube 0.5 mL</p> </div> </div>
Sample Handling	<p>Gently invert tube 5-10X after collection.</p> <p>Let sample clot for approximately 15-20 minutes at room temperature.</p> <p>Centrifuge @ approximately 2500 rpm for 10-15 minutes.</p>	<p>Fill tube to appropriate fill line. (Proper blood to anticoagulant ratio is important for accurate results.)</p> <p>Gently invert the tube 5-10X to prevent clotting.</p>
Sample Storage <i>Note: analyte stability depends on storage temperature and biological half-life.</i>	<p>Draw off serum into a transport tube for storage and shipment. Label the tube with sample ID and collection date and sample type.</p> <p>Store samples refrigerated @ 2-8°C if shipping within 24 hours of collection.</p> <p>Otherwise, freeze samples in transport tubes at -20°C .</p>	<p>Refrigerate. Do Not Freeze.</p>
Sample Shipping	<p>Ship samples on cold packs or dry ice.</p>	<p>Ship samples on cold packs.</p> <p>Make sure tubes do not come in direct contact with cold packs during shipment.</p>

*Do not use SST for progesterone or therapeutic drug monitoring. Use RTT with clot activator instead.

IDEXX BioAnalytics provides these tubes at no charge. To request tubes or for any other questions contact Client Support Services at 1-800-669-0825 or idxxbioanalytics@idexx.com

Clinical Pathology Tube Guide

IDEXX BioAnalytics provides these tubes at no charge. To request tubes or for any other questions contact the Research Support Team at 1-800-669-0825 or idxxbioanalytics@idexx.com

Sample Type	Citrated Plasma Coagulation (ex. PT, PTT, fibrinogen)	Sample Transfer Tubes	Urinalysis	Microbiology
Tube Types	 <p>BTT- Blue Top Tube</p> <p>1.8 mL or 2.7 mL</p> <p>Light blue top with a black ring. sodium citrate anticoagulant (3.2%).</p>	 <p>NON-ADDITIVE FOR TRANSPORT Tube</p> <p>1.5 mL or 6 mL</p> <p>White cap with a black ring. Plain, silicone-coated plastic tube.</p>	 <p>NON-ADDITIVE Tube</p> <p>1.5 mL or 6 mL</p> <p>White cap with a black ring. Plain, silicone-coated plastic tube.</p>	<p>Refer to the Microbiology Sample Preparation and Handling Guide</p>
Sample Handling	<p>Fill tube to appropriate volume (proper blood to anticoagulant ratio is important for accurate results). Gently invert tube 5-10X to prevent clotting. Centrifuge as soon as possible @ 1500 rpm for 15 minutes. Transfer the plasma into a sample transfer tube. Label the tube with sample ID, collection date and "citrated plasma."</p>	<p>Label the tube with sample ID, collection date and sample type (ie. "serum").</p>		
Sample Storage <i>Note: analyte stability depends on storage temperature and biological half-life.</i>	<p>Freeze at -20°C until shipped for analysis.</p>	<p>Refer to the sample handling, storage and shipment directions for the appropriate sample type.</p>	<p>Refrigerate, do not freeze.</p>	
Sample Shipping	<p>Ship samples on cold packs or dry ice.</p>		<p>Ship on cold packs.</p>	

IDEXX BioAnalytics provides these tubes at no charge. To request tubes or for any other questions contact Client Support Services at 1-800-669-0825 or idxxbioanalytics@idexx.com